



Simply the best!

This non-sterile glove is made with a strong durability to provide excellent performance to protect clinicians from bacteria, viruses, and chemicals during procedures where the risk of fluid exposure is moderate to high.

KEY FEATURES AND BENEFITS:

- Ideal for moderate to high risk settings.
- Highly durable synthetic, eliminating the problem of exposure to latex exam gloves while maximizing protection.



Cobalt Blue



Light Purple



Medical



Laboratory



Industry



Other

- LOW DERMATITIS
- AMBIDEXTROUS
- HIGH QUALITY

PRODUCT FEATURES

- The product is for single use only.
- Does not contain natural rubber latex.
- Odor-free and powder-free.
- Can withstand IPA washing (ISO Propyl Alcohol).
- Can be used in medical, dental, patient care, etc.
- Ideal for professional health workers.
- ASTM D-6978



SHIELDGUARD™

PROFESSIONAL GLOVES SERIES

Extraordinary premium gloves with trusted and advanced industry standards and quality

Medical / ASTM D-6978 Series

PREMIUM 550 MEDICAL NITRILE GLOVES

FEATURES

Material	Nitrile
Type	Medical / Chemotherapy Tested
Sterile	No
Powder Free	Yes
Latex Free	Yes
Beaded Cuff	Yes
Ambidextrous	Yes
Color	Cobalt Blue / Light Purple
Outer Surface	Finger Textured
Length*	242 mm
Weight (M)*	4.2g ± 0.2
Single Wall	Yes
Thickness - Palm*	4.7 mil+/- 1.0
Thickness - Fingers*	5.0 mil+/- 1.0
Sizes	S / M / L / XL

STANDARDS & ACCREDITATIONS

FDA 510K
ASTM D6319 / ASTM D6978
EN 374, EN 455
CE
ISO 9001, ISO 13485
GMP

PACKAGING

Box Size	100 Pcs / Box
Case Size	10 Boxes / Case(0.58 kg)

MANUFACTURED

Manufactured for	Shieldguard, Ltd.
Country of Origin	Thailand, Malaysia

*The length, weight and the thickness can be slightly vary depending on the lot numbers.



SHIELDGUARD™ CHEMOTHERAPY GLOVES

PREMIUM 550 MEDICAL NITRILE GLOVES

- LOW DERMATITIS
- HIGH TENSILE FORCE
- ASTM D-6978 TESTED

Chemotherapy

Designed to reflect the demands of the most exacting of industry sectors and benefiting from SHIELDGUARD's proprietary extended hand-created manufacturing process, this range of gloves are significantly easier to don and remove without irritating the skin.

SHIELDGUARD's product line has been proven for over 20 years by many of the professional industry leaders. Only a glove for exposure with superior quality, tensile strength, and best fit. SHIELDGUARD is recognized as the most consistently preferred and best quality glove on the market.



WE PROTECT PEOPLE
SHIELDGUARD

PREMIUM 550 MEDICAL NITRILE GLOVES

- NON-STERILE
- POWDER-FREE
- ASTM D-6978 TESTED

100 Gloves
90 Gloves XL

WE PROTECT PEOPLE
SHIELDGUARD



- LOW DERMATITIS
- AMBIDEXTROUS
- HIGH TENSILE FORCE
- ASTM D-6978 TESTED

PREMIUM 550 MEDICAL NITRILE GLOVES

Powder-Free Non-Sterile Gloves



100 Gloves S-L
90 Gloves XL

PREMIUM 550 MEDICAL NITRILE GLOVES

- NON-STERILE
- POWDER-FREE
- ASTM D-6978 TESTED

100 Gloves
90 Gloves XL



Extraordinary premium gloves with trusted and advanced industry standards and quality

PREMIUM 550 MEDICAL NITRILE GLOVES

- LOW DERMATITIS
- HIGH TENSILE FORCE
- ASTM D-6978 TESTED

Chemotherapy

- S
- M
- L
- XL



WE PROTECT PEOPLE
SHIELDGUARD™
STRENGTH, DURABILITY, COMFORT & FIT

PREMIUM MEDICAL NITRILE GLOVES

- NON-STERILE
- POWDER-FREE
- SINGLE USE ONLY

PRODUCT FEATURES

- The product is for single use only.
- Does not contain natural rubber latex.
- Odor-free and powder-free.
- Can withstand IPA washing ISO Propyl Alcohol.
- Can be used in food industry.
- Ideal for professional health workers.
- ASTM D-6978 Tested

DIRECTION

Dry hands before donning the gloves. Take out the product after checking package. Gloves need to be donned carefully, do not apply excessive force to pull the gloves when donning. Keep the fingertips short and change gloves frequently to avoid bearing issue. The gloves should be discarded with an appropriate waste disposal container after use.

STORAGE CONDITION

The atmosphere must always be conducive to the storage of the gloves. Gloves before and after opened for use shall be stored in dark, cool and dry place. Gloves need to avoid area with excessive heat. The storage areas need to be controlled below 40°C / 104°F. Gloves need to be shielded from direct sunlight and fluorescent lighting.



PRODUCT OF
ShieldGuard
7440 Santa Monica Blvd., Suite 301
Beverly Hills, CA 90210 USA
www.ShieldGuardGlove.com

DISTRIBUTED BY
Shield & Protect, Ltd.
16 Hordle Gardens, At Albans
England AL1 1JW
www.shield_protect.co.uk

Caution: 1. Stay away from sharp objects to avoid puncture. 2. When done using, do not touch the outside of the gloves and dispose of them in a trash can.
Allergy Information: Gloves may cause allergic reactions for some individuals. Do not expose these gloves to individuals with a known or suspected sensitivity to any components used in manufacturing these gloves before consulting a physician. If allergic reactions occur, discontinue use immediately and consult with a physician.

MADE IN THAILAND

LOT Lot/batch number

Date of manufacture

Expiry date

Valid for 3 years from the date of manufacturing



22 X 12 X 7cm

Estimate Size: Large
100 Gloves Per Dispenser Box
10 Boxes/ Carton
Box Dimension: 22 X 11 X 7cm
Carton Dimension: 35.5 X 24.5 X 22.5cm
Carton Weight: 5.6kg
3,300 Cartons/ 40ft He 63CBM
1 Pallet (W: 122cm / D: 102cm) = 72 Cartons
Weight (Full Pallet): 398kg



WE PROTECT PEOPLE
SHIELDGUARD
NITRILE GLOVES NON-STERILE

Inquiries: www.ShieldGuardGlove.com

Product by **McKINLEY**